

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)

Sample : Anterior nasal, Nasal, Throat, saliva swab



Product Features

- High Accuracy, Specificity and Sensitivity
- No need instrument, get results in 15 minutes
- Room temperature storage
- Sample : Anterior nasal, Nasal, Throat, Saliva
- Detect the presence of viral proteins
- Identify acute or early infection

Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold)


Sample Collection

- 1** Please wash and dry your hands thoroughly before the test.


- 2** Read the Instructions for Use carefully.


- 3** Please scan the QR code to watch the operation video.


- 4** Take the swab out of the package and do not touch the sampling end.


- 5 Sample collection:**

 - (1) Anterior nasal swab: Carefully insert the swab 1.5 cm into the nostril until resistance is felt. Under moderate pressure, swab the surface of the nostril 4 to 6 times for at least 15 seconds.




Repeat sampling with the same swab in the other nostril.


 - (2) Nasal swab: The sampler holds the swab to enter the nostril, and when the tip of the swab reaches the back wall of the nasopharyngeal cavity, gently rotate it for a circle, and then slowly take out the swab.


 - (3) Throat Swab: The swab crosses the base of the tongue, and wipes the tonsils back and forth with slight force on both sides of the person being collected for at least 3 times, and then wipes up and down the posterior pharyngeal wall at least 3 times.


 - (4) Saliva swab: Please refrain from eating any food or drink, including chewing gum or tobacco products, within 30 minutes prior to sampling. The saliva specimen is collected by pressing the tip of the tongue against the base of the mandible to collect saliva. Place the sterile disposable swab under the tongue for at least 40 seconds and rotate it 5 or more times to soak completely. Note: If saliva is not collected properly, false negative results may occur.


- 6** Open the larger end of the sampling tube, and insert the swab after collecting the sample into the sampling tube.



- 7** Soak the swab in the buffer for at least 15 seconds, stir the swab several times, and squeeze the swab head 3 times.


- 8** When remove the swab, squeeze the sides of the sampling tube.


- 9** Close the sampling tube with the tube lid.


- 10** Open the foil pouch and place the test cassette on a flat surface.


- 11** Open the tube lid at the front end of the sampling tube, and drop 4 drops of the treated sample into the sample well (S) of the test cassette.



- 12** After reacting at room temperature for 15 minutes, observe the results. A result after 30 minutes is invalid.


- 13** After the test, pack all the components of this test in the Biohazard specimen bag for contaminated waste and dispose of this bag closed with the residual waste. Not reusable.


- 14** Wash or disinfect hands again.



Interpretation of result



Product information

| Product name | Test samples | Specifications | Storage conditions |
|---|--------------------------------------|----------------------------------|--------------------|
| Novel Coronavirus 2019-nCoV Antigen Test (Colloidal Gold) | Anterior nasal, Nasal, Throat Saliva | 1T/kit, 5T/kit, 20T/kit, 25T/kit | 4-30°C |