

- Monkeypox virus antigen can be detected simultaneously.
- No need instrument, get results in 15 minutes.
- Room temperature storage.
- Multiple sample types can be used for detection

# Monkeypox Virus Antigen Rapid Test

## Sample types:

human serum, plasma, whole blood, throat swab, saliva,  
lesion exudate or scab samples

## Model:

Model A: 1T/kit, 5T/kit, 20T/kit, 25T/kit, 40T/kit, 50T/kit

Model B: 1T/kit, 5T/kit, 20T/kit, 25T/kit, 40T/kit, 50T/kit

